



Drug release tests

UNIVERSITY OF PÉCS
INSTITUTE OF PHARMACEUTICAL TECHNOLOGY AND
BIOPHARMACY

LADME system





1897

- Noyes & Whitney

- "The Rate of Solution of Solid Substances in Their Own Solution."
- Rate of drug dissolution is regulated by the saturated solution on the surface of the API particle

1900

- Brunner & Tolloczko

- Prove that the dissolution depends on the physico-chemical character of the API particle (surface, temperature...)



1904

- Nernst & Brunner

- Modify the Noyes-Whitney formula by integrating Fick's law of diffusion

IVIVC

- FDA – in vivo-in vitro correlation – mathematical aspect

1934

- **Pharmacopoea Helvetica** – describes the disintegration test for the first time




1970

- USP 18 introduced the first official drug dissolution test: rotating basket method (Apparatus 1)

2003

Volume 10, Issue 1 | FEBRUARY 2003



Dissolution Technologies

6 FIP/AAPS Guidelines for Dissolution/~~In Vitro Release~~ Testing of Novel/Special Dosage Forms

Martin Siewert,¹ Jennifer Dressman,² Cynthia Brown,³ Vinod Shah⁴

¹ Aventis, Frankfurt, Germany

² JW Goethe University, Frankfurt, Germany


³ Eli Lilly and Company, Indianapolis, IN

⁴ Office of Pharmaceutical Science, Center for Drug Evaluation and Research,
Food and Drug Administration, Rockville, MD

Table 1: Apparatus used for Novel/Special Dosage Forms

| Type of Dosage Form | Release Method |
|--|---|
| Solid Oral Dosage Forms (conventional) | Basket, Paddle, Reciprocating Cylinder or Flow Through Cell |
| Oral Suspensions | Paddle |
| Oral disintegrating Tablets | Paddle |
| Chewable Tablets | Basket, Paddle or Reciprocating Cylinder with glass beads |
| Transdermals – Patches | Paddle Over Disk |
| Topicals – Semisolids | Franz Cell Diffusion System |
| Suppositories | Paddle, modified Basket or Dual Chamber Flow Through Cell |
| Chewing Gum | Special apparatus (Ph.Eur.) |
| Powders and Granules | Flow Through Cell (powder/granule sample cell) |
| Microparticulate Formulations | Modified Flow Through Cell |
| Implants | Modified Flow Through Cell |


Disintegration tester





Disintegration tester

Test for suppositories



Drug release of dosage forms

- During drug release the appearance of the API is determined in the dissolution media in function of time.

According to the equipment structure:

- Closed system
- Open system
- Diffusion system

According to the changing amount of API in the dissolution media:

- Cumulative tests
- Differential tests

Conditions of the test

- Applied equipment
- Dissolution medium
 - Distilled water
 - HCl
 - Buffer
 - Artificial gastric juice
 - Artificial intestinal juice
 - Surfactants
- Rotation per minute
- Sampling time, method and amount
- Analytics

Dissolution tests

Recommended pH values

| pH | Dissolution media |
|----------------|-----------------------------|
| pH 1.0 | HCl |
| pH 1.2 | HCl & NaCl |
| pH 1.5 | HCl & NaCl |
| pH 4.5 | Phosphate or acetate buffer |
| pH 5.5 and 5.8 | Phosphate or acetate buffer |
| pH 6.8 | Phosphate buffer |
| pH 7.2 and 7.5 | Phosphate buffer |



| 0-1 h | 1-2 h | 2-3 h | 3-4 h | 4-5 h | 5-6 h | 6-7 h | 7 h |
|-------|-------|-------|-------|-------|-------|-------|-----|
| 1.0 | | | | | | | |
| 1.2 | 6.8 | | | | | | |
| 1.2 | 2.5 | 4.5 | 7.0 | | 7.5 | | |
| 1.5 | 4.5 | | | 7.2 | | | |

Dissolution test

Increasing the IVIVC correlation


| Physiological factors | In vitro factors |
|------------------------------------|----------------------------------|
| pH | Different dissolution media (pH) |
| GI motility | Stirring conditions |
| Fat and protein, food interactions | Adding fat or milk |
| Enzymes | Applying enzymes |
| Bile | Applying surfactants |
| GI transit time | pH-gradient tests |

Dissolution testers (1,2,3,4)




Dissolution test

Solid dosage forms – rotating basket apparatus



Dissolution test

Solid dosage forms – paddle apparatus





International Conference on Harmonisation (ICH)

- Rotating basket method (50/100 rpm)
- Rotating paddle method (50/75 rpm)
- Sampling time every 15 min, at rapid preparations in every 5-10 min
- Volume: 500 ml, 900 ml, 1000 ml (sink conditions)
- pH=1.2-6.8 (not more than pH=8.0)

Dissolution run


Required Steps



Requires user presence and interaction (technician-dependent)

Does not require user presence and interaction (technician-independent)

Drug dissolution procedure



Degassing

- Incorrect degassing can lead to decreased or increased API levels:
 - Bubbles can adhere on the surface of the preparation (tablet) and it can swim on the dissolution medium surface
→ different hydrodynamic conditions
 - Bubbles can adhere on the surface of the rotating basket

Dissolution tester



Dissolution tester



Sampling canulla



UHMW – poly-ethylene – 35-40 μ m average pore size

Canulla filters



Paddle



PTFE – paddle



Basket



Basket



Basket






Sinker




Sinker



Place of sampling


- It should be taken from the upper level
- It should be taken from the middle level
- Sampled at the bottom level
 - injected
 - substituted
 - taken in




from level and
container.
be

Flow-through cell

- Structure:
 - **Container** of the medium
 - **Pump** – pumps the medium through the cell
 - **Flow-through cell**
 - **Water bath**




Transdermal drug dissolution







Franz-cell



MicroettePlus






Transdermal patches dissolution test

- Disk assembly method
 - Rotating paddle is used
 - 125 μm steel sieve is applied
 - Patch should be placed in the container with its release surface facing up
 - Temperature: $32 \pm 0.5^\circ\text{C}$

Dissolution tests for patches

Disk assembly method






USP Apparatus 5

Transdermal patches dissolution test

- Cell method
 - Paddle is used
 - Cell is stainless steel

Dissolution tests for patches

Cell method




Transdermal patches dissolution test

- Rotating cylinder method
 - Paddle is substituted with a cylinder

Dissolution tests for patches

Rotating cylinder method




USP Apparatus 6

Chewing gum's dissolution test


- Buffer solution medium- artificial saliva (20-40 ml), is placed into a chewing chamber
- Temperature: $37 \pm 0.5^{\circ}\text{C}$
- Chewing frequency: 60 cycles/min

Dissolution test for chewing gums





Kinetic models








Evaluation of drug release tests

UNIVERSITY OF PÉCS
INSTITUTE OF PHARMACEUTICAL TECHNOLOGY AND
BIOPHARMACY

Liberation

Types of drug release curves



Evaluation of drug release

- 1.) model-independent
- 2.) model-dependent evaluation

Evaluation of drug release

Model-independent method

- Fitted functions sometimes have no theoretical relationship with dissolution, they are based on experiences
- Independent from kinetic models

Evaluation of drug release

Model-independent method

Used to evaluate and for the comparison of different drug dissolution profiles

Two types:

1. Characterization of specific values of the dissolution profiles (mean dissolution time, area under the curve, dissolution efficiency) – statistical evaluation/comparison


Model-independent methods

Dissolution Efficiency (*D.E.*)

Ratio of AUC of dissolution curve and the AUC of the 100% dissolution.

$$D.E. = \frac{\int_0^t m_t dt}{m_{100\%} t} \times 100$$

m_t cumulative dissolution at ' t ' time in %
 $m_{100\%}$ 100% dissolution



Model-independent methods

Mean Dissolution Time (MDT)

Time of the 63,2% dissolution.

$$MDT = \frac{\int_{i=1}^n \bar{t}_i \Delta m_i}{\int_{i=1}^n \Delta m_i}$$

$$\bar{t}_i = \frac{t_{i-1} + t_i}{2}$$

- i no. of samples,
- n no. of sampling times,
- \bar{t}_i average time between t_i and t_{i-1}
- Δm dissolved API amount between t_i and t_{i-1}

Evaluation of dissolution tests

Model-independent method

2. Pairwise comparison of dissolution data by appropriate statistical methods

Model-independent methods

fit factors (f_1, f_2)

Method compares the reference dissolution data to the test preparation .

f_1 difference factor

$$f_1 = \frac{\sum_{i=1}^n (R_i - T_i)}{\sum_{i=1}^n R_i} \times 100$$

f_2 similarity factor

$$f_2 = 50 \lg \left\{ \left[1 + \frac{1}{n} \sum_{i=1}^n (R_i - T_i) \right]^{-0,5} * 100 \right\}$$

n no. of samples,
 R_i dissolution % at 'i' time of the reference preparation,
 T_i dissolution % at 'i' time of the test preparation

Model-independent methods

fit factors (f_1, f_2)

If the *diference* (f_1) factor =0, then the dissolution of the reference and the test preparation are considered to be equal.



Model-independent methods

fit factors (f_1, f_2)

Similarity (f_2) factor is a figure between 0 and 100.
If $f_2 > 50$, then the two curves are similar.

Evaluation of drug release tests

Model-dependent method


1.) With theoretically valid background
2.) Without theoretical background, based on experiences

Model-dependent methods

Zero order dissolution kinetics

$$m = kt$$

$$c = \frac{m}{V}$$



m API amount at 't' time

t time

c API % at 't' time

k rate constant


Model-dependent methods

First order kinetics

$$\frac{dm}{dt} = km$$

$$m = m_0 e^{-kt}$$

$$(\ln m - \ln m_0) = -kt$$




Model-dependent methods

polimer matrix diffusion

Higuchi formula

$$m = \sqrt{D(2C - C_s)C_s t}$$

m *dissolved API amount at 't' time*
 D *diffusion constant*
 C_s *API solubility*
 C *API concentration at zero time*




Model-dependent methods

porous systems

$$m = \sqrt{\frac{D_s \varepsilon C_s}{\tau} (2M - \varepsilon C_s) t}$$

- m *dissolved API % at 't' time*
 M *API amount in the matrix*
 C_s *API solubility*
 D_s *API diffusion constant*
 ε *porosity*
 τ *convolution constant*



Model-dependent methods

Biological erosion

$$m = k_e A_e t$$

m *dissolved API % at 't' time*
 k_e *erosion constant*
 A_e *eroding (changing) surface*

Model-dependent methods

Osmotic systems

$$m = \frac{\pi A}{h} \Delta p C_s t$$

m *dissolved API % at 't' time*

π *permeability*

A *surface*

$p\Delta$ *difference in osmotic pressure*

C_s *solubility*

h *thickness of the semipermeable membrane*

Model-dependent methods

General model: Weibull model

$$m / m_{\infty} = 1 - \exp\{-[(t - t_0) / \tau]^{\beta}\}$$

$$m / m_{\infty} = 1 - e^{-\left[\frac{t-t_0}{\tau}\right]^{\beta}}$$

m dissolved API % at 't' time


m_{∞} dissolved API % at $t = \infty$ time

t time

t_0 lag time

τ time of 63,2% API dissolution


β shape factor



Model-dependent methods

Weibull model

β shape factor



$\beta=1$ – ZOK dissolution

Solid oral dosage forms

Immediate release typically means that 75% of the API is dissolved within 45 minutes

- Rapidly dissolving: $\geq 85\%$ in ≤ 30 minutes
- Very rapidly dissolving: $\geq 85\%$ in ≤ 15 minutes

Comparative dissolution testing

Profile similarity determination

Two conditions to determine if the dissolution profiles of two products/batches in a particular dissolution medium are similar:

1. If **both the test and reference product show more than 85% dissolution within 15 minutes**, the profiles are considered to be **similar**
 - **No calculations are required**

If this is not the case, apply point 2

2. Calculate the f_2 value (similarity factor):
 - If **$f_2 \geq 50$** , the profiles are **normally regarded similar**

Comparative dissolution testing

Similarity factor f_2

$$f_2 = 50 \cdot \log \left[\frac{100}{\sqrt{1 + \frac{\sum_{t=1}^{t=n} [\bar{R}(t) - \bar{T}(t)]^2}{n}}} \right]$$

n = number of time points

$R(t)$ = mean % API dissolved of **reference product** at time point x

$T(t)$ = mean % API dissolved of **test product** at time point x

- **Minimum of 3 time points** (zero excluded)
- **12 units (each in own dissolution vessel) for each product** (for “official” purposes)
- **Only one measurement should be considered after both products have reached 85 % dissolution**
- **RSD** at higher time points $\leq 10\%$

Comparative dissolution testing

Dissolution conditions (**study design**)

| | |
|--|---|
| Apparatus (choice) | <ul style="list-style-type: none"> • Paddle, 50 (75) rpm or • Basket, 100 rpm |
| Dissolution media All three media for full comparison | <ol style="list-style-type: none"> 1. Buffer pH 6.8 <u>or</u> simulated intestinal fluid without enzymes 2. Buffer pH 4.5 3. 0.1 M HCl <u>or</u> buffer pH 1.2 <u>or</u> simulated gastric fluid without enzymes |
| Volume of media | 900 ml or less |
| Temperature | 37°C ± 0.5°C |
| Sampling points | 10, 15, 20, 30, 45, (60, 120) min. (typical) |
| Units (individual) | 12 for “official” studies |

Typical time points

Immediate release tablets (capsules)

| Point | Time |
|-------|------|
| 1 | 10 |
| 2 | 15 |
| 3 | 20 |
| 4 | 30 |
| 5 | 45 |

Rationale:

1. Condition 1

- $\geq 85\%$ dissolution of both products within 15 minutes
- 15 minute time point thus **essential**

2. Condition 2, for calculation of f_2

- a **minimum of 3 points** are required
- Only **one measurement should be considered after 85 % dissolution** (both tablets)
- **20 minute time point thus first possible one** (if 15 minute fails 1st condition)

Example

Determination of similarity of profiles

| Example 1-A | | |
|----------------|---------------------------|-----------------|
| Time (min) | % API dissolved | |
| | Tablet A (Ref) | Tablet B (Test) |
| 10 | 87 | 94 |
| 15 | 96 | 99 |
| 20 | 99 | 99 |
| 30 | 100 | 99 |
| 45 | 101 | 99 |
| 60 | 101 | 99 |
| f2 required? | No, $\geq 85\%$ in 15 min | |
| f2 (n = N/A ?) | profiles similar | |

| Example 1-B | | |
|--------------|-----------------|-----------------|
| Time (min) | % API dissolved | |
| | Tablet D (Ref) | Tablet E (Test) |
| 10 | 55 | 57 |
| 15 | 72 | 78 |
| 20 | 85 | 91 |
| 30 | 97 | 100 |
| 45 | 102 | 100 |
| 60 | 103 | 101 |
| f2 required? | Yes | |
| f2 (n = 3 ?) | 64 (similar) | |

Example

Determination of similarity of profiles (cont.)

| Example 1-C | | |
|--------------|-----------------|-----------------|
| Time (min) | % API dissolved | |
| | Tablet X (Ref) | Tablet Y (Test) |
| 10 | 29 | 34 |
| 15 | 38 | 41 |
| 20 | 47 | 50 |
| 30 | 63 | 64 |
| 45 | 80 | 79 |
| 60 | 95 | 91 |
| f2 required? | Yes | |
| f2 (n = 6 ?) | 74 (similar) | |

| Example 1-D | | |
|--------------|------------------|-----------------|
| Time (min) | % API dissolved | |
| | Tablet A (Ref) | Tablet Y (Test) |
| 10 | 87 | 55 |
| 15 | 96 | 72 |
| 20 | 99 | 85 |
| 30 | 100 | 97 |
| 45 | 101 | 102 |
| 60 | 101 | 103 |
| f2 required? | Yes | |
| f2 (n = 3 ?) | 31 (not similar) | |

Thank you for the attention!